

# TUMOR MARKERS



## Diagnosis & Monitoring of **CANCER**

from diagnosis,  
the seeds of better health



# TUMOR MARKERS

## TUMOR MARKERS


### ■ From dreams...

An ideal tumor maker (**TM**) would be a substance produced only by the neoplastic cells, specific to one type of tumor (no false positives) and detectable right from the initial stage of the disease (no false negatives).

It would be undetectable in healthy subjects, and enable the screening and diagnosis of cancer. The [TM] level would correlate closely with tumor size, contribute to the initial extension of the profile and evaluation of therapeutic efficacy, as well as the early detection of recurrent diseases.

### ■ ...to reality

Most TMs are also produced by normal cells, and have only average sensitivity and specificity levels which depend on the decision-making threshold chosen.




This booklet was created with the collaboration of :

Doctor RIEDINGER

Laboratoire de biologie,

Centre G.F. Leclerc-CRLCC, Dijon, FRANCE

The TM level is influenced by a large number of parameters.



Evolution of TM level over time is closely correlated to variations in the size of the tumor.

### Non-specific increases in TM levels other than tumor pathologies

<b>CEA</b>	Smoking, benign pulmonary and digestive pathologies, chronic renal insufficiency
<b>CA 19-9</b>	Benign pulmonary and digestive pathologies (cholecystitis and pancreatitis ++)
<b>AFP</b>	Benign hepatic pathologies, multiple pregnancies
<b>PSA</b>	Prostate adenoma, prostatitis, prostate manipulations
<b>βhCG</b>	Pregnancy
<b>CA 15-3</b>	Benign hepatic, digestive, pulmonary and ovarian pathologies, auto-immune disease
<b>CA 125</b>	Serosa pathologies, cirrhosis, ascitis, endometriosis, end of pregnancy
<b>CT</b>	Apudoma, chronic renal insufficiency, hyperparathyroidism, Paget's disease
<b>Tg</b>	Benign thyroid diseases
<b>NSE</b>	Apudoma, renal tumors, benign pulmonary pathologies
<b>β2-M</b>	Rheumatoid arthritis, AIDS, lupus, crohn's disease, renal tubular dysfunction.

## CLINICAL USE OF MARKERS

Organ	Histological type	Main or secondary marker
<b>Breast</b>	Adenocarcinoma	CA 15-3 - CEA
<b>Ovary</b>	Serosa	CA 125
	Mucosa	CA 19-9 (CEA - CA 72-4)
	Germinal	AFP - βhCG
	Granulosa	Inhibin (estradiol)
<b>Uterus (neck)</b>	Epidermoid	SCC or CYFRA 21-1
<b>Endometrium</b>	Adenocarcinoma	CA 125 - CA 19-9
<b>Placenta</b>	Trophoblastic	βhCG
<b>Prostate</b>	Adenocarcinoma	PSA - free PSA (PAP)
<b>Bladder</b>	Adenocarcinoma	TPA
<b>Testis</b>	Seminoma	βhCG
	Non seminoma	βhCG - AFP
<b>Oesophagus</b>	Epidermoid	SCC or CYFRA 21-1
	Adenocarcinoma	CA 19-9 (CEA)
<b>Stomach</b>	Adenocarcinoma	CA 72-4 (CA 19-9 - CEA)
<b>Liver</b>	Hepatocarcinoma	AFP - DCP (CEA)
	Metastasis (unknown primitive)	CEA - AFP - CA 19-9 - CA 15-3 - NSE
<b>Bile ducts</b>	Adenocarcinoma	CA 19-9 (CEA)
<b>Pancreas</b>	Adenocarcinoma	CA 19-9 (CEA - CA 50)
	Endocrine	Digestive hormones - NSE
<b>Small intestine</b>	Carcinoid	5-HIAA* (CEA - CA 19-9)
<b>Colon/rectum</b>	Adenocarcinoma	CEA (CA 19-9)
<b>Anus</b>	Epidermoid	SCC or CYFRA 21-1
<b>Medulo-adrenal</b>	Pheochromocytoma	VMA* - metanephrines* - catecholamines*
<b>Lung</b>	Adenocarcinoma	CEA
	Epidermoid	SCC or CYFRA 21-1
	Small cells	NSE
<b>ENT</b>	Epidermoid	SCC or CYFRA 21-1
<b>Thyroid</b>	Medullary	Calcitonin - CEA
	Differentiated	Thyroglobulin
<b>Nervous system</b>	Neuroblastoma	NSE - HVA* - catecholamines*
<b>Hemato-lymphoid</b>	Lymphoid malignancies	β2-M

\* assays using urine

## The clinical use of TMs has 5 main indications

CEA	.....	Carcinoembryonic Antigen
AFP	.....	$\alpha$ -Fetoprotein
CA	.....	Cancer Antigen
$\beta$ 2-M	.....	$\beta$ 2-Microglobulin
CT	.....	Calcitonin
5-HIAA	.....	5-Hydroxyindoleacetic Acid
5'Nu	.....	5'Nucleotidase
DCP	.....	Decarboxyprothrombin
$\gamma$ GT	.....	$\gamma$ -Glutamyl Transferase
hCG	.....	Human Chorionic Gonadotropin
HVA	.....	Homovanillic Acid
LDH	.....	Lactate Dehydrogenase
NSE	.....	Neurone Specific Enolase
ALP	.....	Alkaline Phosphatase
BAP	.....	Bone Alkaline Phosphatase
PAP	.....	Prostatic Phosphatase Acid
PSA	.....	Prostate Specific Antigen
SCC	.....	Squamous Cell Carcinoma
Tg	.....	Thyroglobulin
TPA	.....	Tissue Polypeptide Antigen
VMA	.....	Vanillylmandelic Acid

## SCREENING

PSA is the only maker that is used for screening of prostatic cancer in conjunction with digital rectal examination in males above the age of 50. In addition, **targeted screening** of thyroid medullary cancer using CT and hepatocarcinoma using AFP is justified.

## DIAGNOSIS

The diagnostic use of TM is rare and is limited to the following assays :

- $\beta$ hCG for germinal tumors,
- AFP for hepatocarcinomas,
- CT for thyroid medullary cancer,
- PSA for prostate cancer.■

■ *To improve the diagnostic performance of PSA, different criteria have been proposed :*  
*The free PSA/total PSA ratio,*  
*the PSA density (ratio of PSA to gland size),*  
*the PSA velocity (annual increase of PSA)*  
*and the level of PSA linked to patient age.*

## PROGNOSIS

The prognostic value of TMs is generally linked to the tumor size and often indicates the extent of disease. The only TMs with an independant prognostic value for the extent of disease are :

- CEA in colorectal cancer,
- $\beta$ hCG and AFP for germinal tumors,
- CYFRA 21-1 in non-small cell pulmonary cancer,
- LDH in germinal tumors and lymphoma.

## THERAPEUTIC MONITORING

Therapeutic monitoring is evaluated by tracing the evolution of TM levels measured using the same technique within the same laboratory. Each patient is considered to be their own control and any reference to normal values should be abandoned.

### Individual kinetics of the marker level

Integration of clinical and therapeutic data to the TM individual kinetics graph, as logarithmic coordinates, enables the different kinetic parameters to be calculated.

To evaluate therapeutic efficacy, the TM half-life ( $T_{1/2}$ ) and its minimum level obtained during treatment (min[TM]) are measured.

### Apparent half-life ( $T_{1/2}$ )

The half-life is the time taken for the TM concentration to decrease by half in the initial phase of therapy. This parameter measures the rate at which the size of the tumor size is reduced.

Its value depends on the type of therapy administered, as well as its efficacy. After radical surgery, the  $T_{1/2}$  reflects only the clearance of TM.

While administering chemotherapy, radiotherapy or hormone therapy, the  $T_{1/2}$  also reflects the transitory residual tumor secretion.

- *Some paradoxical increases in the TM level in the initial phase of efficient systemic therapy result from massive cytolysis or regeneration of hepatocytes.*

*When a heterogenous tumor secretes several TMs, a discrepant evolution of these markers can indicate a different sensitivity of the secreting cell clones to the therapy administered.*

## Minimum level obtained during therapy min[TM]

The min[TM] is a good indicator of residual disease. Its value and the time taken to obtain the result depends on the type of therapy administered and its efficacy.

- *The persistence of a high min[TM] level reveals the existence of a remaining secreting tumor. A return to usual values is not always associated with sterilization of the tumor. The value can correspond to the persistence of small remaining tumors or the disappearance of the only secreting cell clone.*

*For TMs for which the half-life is short, the persistence of a detectable level 3 weeks after complete surgical resection indicates the persistence of residue, probably due to a tumor.*

## Characteristics of main tumors

TM	$T_{1/2}$ (therapy)	Normal values
CEA	several weeks	< 2.5-5 µg/l
CA 19-9	8.5 days (surgery)	< 37 kU/l
AFP	4-5 days (surgery)	< 10 µg/l
PSA	2.2-3.2 days (surgery) 1.4-2.6 months (radiotherapy)	< 2.5 to 6.5 µg/l*
βhCG	24 to 36 hours (surgery)	< 5-10 U/l
CA 15-3	8-15 days (surgery)	< 25-35 kU/l
CA 125	4.8 days (surgery) 9.2 days (chemotherapy)	< 35 kU/l
CT	20 minutes (surgery)	< 10 ng/l
Tg	2.7 days (surgery)	< 25 µg/l
SCC	20 minutes (surgery)	< 1.5 µg/l
CYFRA 21-1	4 days (surgery)	< 1.8 µg/l
TPA	7 days (surgery)	< 95 U/l

\* 40-49 yrs old < 2.5 µg/l ; 50-59 yrs old < 3.5 µg/l ; 60-69 yrs old < 4.5 µg/l ; 70-79 yrs old < 6.5 µg/l (Standard Hybritech)



# PRACTICAL QUESTIONS

## When to prescribe TM ?

**Before therapy :** to obtain a reference value.

**During therapy :** according to a rhythm adapted to the type of therapy, the half-life and the initial level of TM.

**During therapeutic monitoring :** with a frequency adapted to

- the risk of recurrent disease,
- the average time lapse before the onset of recurrent disease,
- the therapeutic alternatives available.

**During evolving recurrent disease :** according to a rhythm adapted to the doubling time of the TM.

## Analysis of discrepancies

### clinical :

- modification of the filtration functions of TM (hepatic or renal insufficiency),
- tumor heterogeneity,
- massive tumor cytolysis or regeneration of hepatocytes,
- inadequate evaluation (even dissociation) of the clinical response.

### analytical :

- problem of reagent standardization,
- variable sensitivity of reagents to the hook effect or interfering substances,
- problem of international standardization,
- molecular heterogeneity of TM,
- presence of HAMA\*.

\* mouse anti-human immunoglobulin antibodies

## Can TMs be detected elsewhere than in blood ?

After surgery, the expression of TM in tumor cells can be detected in tissue using the immunohistochemical technique.

TMs can also be assayed in biological fluids such as :

- urine (catecholamines and NSE in neuroendocrine tumors),
- pleural fluid (hyaluronic acid in mesotheliomas),
- ascite fluid (CA 19-9 and CEA in peritoneal carcinosis),
- CSF (NSE in organic neurological pathologies),...

## What is the interest of the enzymatic assays in cancerology ?

They can usefully complete the information linked to TMs. If the increase in some enzymatic activities (ie. ALP, 5'Nu,  $\gamma$ GT, LDH) indicates a non-specific response of the liver parenchyma to the invading tumor, other hyperactivities are a more specific result of the existence of a tumor : BAP (primary or metastatic bone tumor), PAP (prostatic tumor), lysozyme (myeloproliferative syndromes) etc...

## ALWAYS REMEMBER

- 1 A normal TM level does not exclude cancer.
- 2 A high TM level does not always indicate cancer.

## ESSENTIAL DO'S...

- 1 Assay the TM level before administering therapy.
- 2 Interpret the TM level taking into account the clinical and radiological context.
- 3 Control all positive results indicating the need for therapeutic decision making using a new sample.
- 4 Monitor patients using the same technique in the same laboratory.
- 5 Store sera in a serum bank.
- 6 Examine the background history when changing assay technique.
- 7 Integrate the individual evolution kinetics of the TM in the patient report.

## ...ESSENTIAL DON'TS

- 1 Assay a marker when no therapy is available.
- 2 Monitor the tumor site using several TMs which are not really complementary.

### Blood Tumor Markers assays :

<b>VIDAS® CEA</b>	ref. 30412
<b>VIDAS® AFP</b>	ref. 30413
<b>VIDAS® TPSA</b>	ref. 30428
<b>VIDAS® FPSA</b>	ref. 30440
<b>VIDAS® CA 15-3®</b>	ref. 30429
<b>VIDAS® CA 19-9™</b>	ref. 30427
<b>VIDAS® CA 125 II™</b>	ref. 30426
<b>VIDAS® β2-M</b>	ref. 30420
<b>VIDIA® TPSA*</b>	ref. 38500
<b>VIDIA® FPSA*</b>	ref. 38510
<b>VIDIA® CEA*</b>	ref. 38520
<b>VIDIA® AFP*</b>	ref. 38530

\* in development

Availability of some VIDAS tests may be restricted in certain countries due to registration requirements. Consult our local representatives for further information.



04-05 / 008GB99067A / This document is not legally binding. bioMérieux reserves the right to modify specifications without notice.  
BIOMÉRIEUX, the blue logo, VIDAS and VIDJA are registered and protected trademarks belonging to bioMérieux sa or one of its subsidiaries  
bioMérieux sa RCS Lyon 673 620 399 / Photos : bioMérieux, Getty images / Printed in France / T.L. McCANN SANTÉ Lyon / RCS Lyon B 398 160 242

**bioMérieux sa**  
69280 Marcy l'Etoile  
France  
Tel. : 33 (0)4 78 87 20 00  
Fax : 33 (0)4 78 87 20 90  
[www.biomerieux.com](http://www.biomerieux.com)

